

Product datasheet for **AM31285FC-N**


MICB Mouse Monoclonal Antibody [Clone ID: B-D54]

Product data:

Product Type:	Primary Antibodies
Clone Name:	B-D54
Applications:	FC
Recommended Dilution:	Flow Cytometry: Use 10 µl to label 10 ⁶ cells or 100 µl of whole blood.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Recombinant human MICB
Specificity:	This antibody recognises both natural and recombinant MICB.
Formulation:	PBS Label: FITC State: Liquid purified IgG fraction Stabilizer: 5% BSA Preservative: 0.09% Sodium Azide Label: Fluoresceinisothiocyanat
Purification:	Ion Exchange Chromatography
Conjugation:	FITC
Storage:	Store the antibody undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.
Gene Name:	MHC class I polypeptide-related sequence B
Database Link:	Entrez Gene 4277 Human Q29980
Synonyms:	MICB, PERB11.2
Protein Families:	Druggable Genome
Protein Pathways:	Natural killer cell mediated cytotoxicity



[View online »](#)

Product images:

A typical staining pattern with the B-D54 monoclonal antibody on Molt4 cell line