

Product datasheet for **AM31286AF-N**

MICB Mouse Monoclonal Antibody [Clone ID: B-D55]

Product data:

Product Type:	Primary Antibodies
Clone Name:	B-D55
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Recombinant human MICB
Specificity:	This antibody recognises both natural and recombinant MICB.
Formulation:	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free State: Azide Free State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Ion exchange chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.
Gene Name:	Homo sapiens MHC class I polypeptide-related sequence B (MICB), transcript variant 1
Database Link:	Entrez Gene 4277 Human Q29980
Synonyms:	MICB, PERB11.2
Protein Families:	Druggable Genome
Protein Pathways:	Natural killer cell mediated cytotoxicity



[View online »](#)