

Product datasheet for **AM31287FC-N**

GNLY Mouse Monoclonal Antibody [Clone ID: B-L38]

Product data:

Product Type:	Primary Antibodies
Clone Name:	B-L38
Applications:	FC
Recommended Dilution:	Flow Cytometry: Use 10 µl to label 10e6 cells or 100 µl of whole blood.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Recombinant human Granulysin (NKG5 and LAG-2)
Specificity:	This antibody recognizes both natural and recombinant Granulysin.
Formulation:	PBS Label: FITC State: Liquid purified Ig fraction Stabilizer: 5% BSA Preservative: 0.09% Sodium Azide Label: Fluorescein isothiocyanat
Purification:	Ion Exchange Chromatography
Conjugation:	FITC
Storage:	Store the antibody undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.
Gene Name:	granulysin
Database Link:	Entrez Gene 10578 Human P22749
Synonyms:	Protein NKG-5, Lymphokine LAG-2, GNLY, LAG2, NKG5, TLA519
Protein Families:	Secreted Protein



[View online »](#)