

Product datasheet for **AM31291AF-N**

Galectin 3 (LGALS3) Mouse Monoclonal Antibody [Clone ID: B-S43]

Product data:

Product Type:	Primary Antibodies
Clone Name:	B-S43
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Recombinant human Galectin-3
Specificity:	This antibody recognises both natural and recombinant Galectin-3 protein.
Formulation:	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free State: Azide Free State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Ion exchange chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.
Gene Name:	lectin, galactoside binding soluble 3
Database Link:	Entrez Gene 3958 Human P17931
Synonyms:	Mac-2, Lgals3, GAL3, GALBP, CBP35, L-31
Protein Families:	Secreted Protein



[View online »](#)