

Product datasheet for **AM31293AF-N**

IL24 Mouse Monoclonal Antibody [Clone ID: B-C60]

Product data:

Product Type:	Primary Antibodies
Clone Name:	B-C60
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Immunogen:	Recombinant human IL-24
Specificity:	This antibody recognises both natural and recombinant IL-24 (interleukin-24) also called MDA-7 or Th2-specific cytokine FISP.
Formulation:	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free State: Azide Free State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Ion exchange chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.
Gene Name:	interleukin 24
Database Link:	Entrez Gene 11009 Human Q13007
Synonyms:	IL-24, MDA7, ST16
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways:	Cytokine-cytokine receptor interaction, Jak-STAT signaling pathway



[View online »](#)