

Product datasheet for **AM31295AF-N**

NCR3 Mouse Monoclonal Antibody [Clone ID: B-S30]

Product data:

Product Type:	Primary Antibodies
Clone Name:	B-S30
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Immunogen:	Recombinant human NKp30/Fc
Specificity:	This antibody recognises both natural and recombinant NKp30.
Formulation:	Phosphate-buffered saline. Sterile filtered through 0.22µm. Carrier and preservative free. State: Azide Free State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Ion exchange chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.
Gene Name:	natural cytotoxicity triggering receptor 3
Database Link:	Entrez Gene 259197 Human O14931
Synonyms:	NKp30, NCR3, 1C7, LY117
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS
Protein Pathways:	Natural killer cell mediated cytotoxicity



[View online »](#)