

## Product datasheet for **AM31296AF-N**

### Mouse IgG1 Isotype Control

#### Product data:

Product Type:	Isotype Controls
Applications:	FC
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Rat Monocytes and Macrophages.
Buffer:	Phosphate-buffered saline. Carrier and preservative free. Sterile-filtered through 0.22 µm and treated to remove endotoxins.
Concentration:	1 mg/ml
Purification:	Ion Exchange Chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C. Shelf life: one year from despatch.



[View online »](#)