

Product datasheet for **AM31296PP-N**

Mouse IgG1 Isotype Control

Product data:

Product Type:	Isotype Controls
Applications:	FC
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Rat Monocytes and Macrophages.
Buffer:	PBS containing 1% BSA and 0.09% Sodium Azide
Purification:	Ion Exchange Chromatography
Conjugation:	PerCP
Storage:	Store the antibody undiluted at 2-8°C. DO NOT FREEZE! Shelf life: one year from despatch.



[View online »](#)