

## Product datasheet for **AM31297BT-N**

### Mouse IgG2a Isotype Control

#### Product data:

Product Type:	Isotype Controls
Applications:	FC
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Rat Monocytes and Macrophages.
Buffer:	PBS containing 5% BSA as stabilizer and 0.09% Sodium Azide as preservative
Reconstitution Method:	Restore with 1.0 ml deionised water.
Purification:	Ion Exchange Chromatography
Conjugation:	Biotin
Storage:	Store the antibody undiluted at 2-8°C after reconstitution. Shelf life: Six months from despatch.



[View online »](#)