

Product datasheet for **AM31297FC-N**

Mouse IgG2a Isotype Control

Product data:

Product Type:	Isotype Controls
Applications:	FC
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Rat Monocytes and Macrophages.
Buffer:	PBS containing 1% BSA and 0.09% Sodium Azide
Purification:	Ion Exchange Chromatography
Conjugation:	FITC
Storage:	Store the antibody undiluted at 2-8°C. Shelf life: one year from despatch.



[View online »](#)