

## Product datasheet for **AM31297PP-N**

### Mouse IgG2a Isotype Control

#### Product data:

Product Type:	Isotype Controls
Applications:	FC
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Rat Monocytes and Macrophages.
Buffer:	PBS containing 1% BSA and 0.09% Sodium Azide
Purification:	Ion Exchange Chromatography
Conjugation:	PerCP
Storage:	Store the antibody undiluted at 2-8°C. <b>DO NOT FREEZE!</b> Shelf life: one year from despatch.



[View online »](#)