

## Product datasheet for **AM31297RP-N**

### Mouse IgG2a Isotype Control

#### Product data:

Product Type:	Isotype Controls
Applications:	FC
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Rat Monocytes and Macrophages.
Buffer:	PBS containing 1% BSA and 0.09% Sodium Azide
Reconstitution Method:	Restore with 1 ml deionised water.
Purification:	Ion Exchange Chromatography
Conjugation:	PE
Storage:	Store the antibody undiluted at 2-8°C after reconstitution. <b>DO NOT FREEZE!</b> Shelf life: Six months from despatch.



[View online »](#)