

Product datasheet for **AM31298FC-N**

Mouse IgG2b Isotype Control

Product data:

| | |
|---------------|---|
| Product Type: | Isotype Controls |
| Applications: | FC |
| Host: | Mouse |
| Isotype: | IgG2a |
| Clonality: | Monoclonal |
| Immunogen: | Rat Monocytes and Macrophages. |
| Buffer: | PBS containing 1% BSA and 0.09% Sodium Azide |
| Purification: | Ion Exchange Chromatography |
| Conjugation: | FITC |
| Storage: | Store the antibody undiluted at 2-8°C. Shelf life: one year from despatch. |



[View online »](#)