

Product datasheet for **AM31384PU-N**

alpha 1 Fetoprotein (AFP) Mouse Monoclonal Antibody [Clone ID: B492M]

Product data:

Product Type:	Primary Antibodies
Clone Name:	B492M
Applications:	ELISA, LF
Recommended Dilution:	ELISA. Lateral flow assays.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Alpha Fetoprotein (AFP)
Specificity:	This antibody reacts to human Alpha Fetoprotein (AFP). Does not crossreact with Human Serum Albumin.
Formulation:	PBS containing 0.09% Sodium azide as preservative State: Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Ion exchange chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for 15 days or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	alpha fetoprotein
Database Link:	Entrez Gene 174 Human P02771
Synonyms:	Alpha-1-fetoprotein, Alpha-fetoglobulin, HPAFP
Protein Families:	Druggable Genome, Secreted Protein



[View online »](#)