

Product datasheet for **AM31388PU-N**

Amphetamine Mouse Monoclonal Antibody [Clone ID: BDI088]

Product data:

Product Type:	Primary Antibodies
Clone Name:	BDI088
Applications:	ELISA
Recommended Dilution:	ELISA.
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Specificity:	This antibody reacts to Amphetamine.
Formulation:	0.05M Phosphate, 0.85% Sodium Chloride, pH 7.2 State: Purified State: Liquid purified Ig fraction Preservative: 0.05% Sodium Azide
Concentration:	lot specific
Purification:	Affinity Chromatography on Protein A
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	Amfetamine



[View online »](#)