

## Product datasheet for **AM31397PU-N**

### **Oxazepam Mouse Monoclonal Antibody [Clone ID: B650M]**

#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	B650M
Applications:	ELISA
Recommended Dilution:	ELISA.
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Benzodiazepine-BSA conjugate
Specificity:	This antibody reacts to Benzodiazepine (BZO).
Formulation:	PBS, pH 7.0 containing 0.02% Sodium azide as preservative State: Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	>95% pure (SDS-PAGE)
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	Benzodiazepine, BZO



[View online »](#)