

## Product datasheet for **AM31429PU-N**

### CTC\_RS14060 Mouse Monoclonal Antibody [Clone ID: B418M]

#### Product data:

|                       |  |
|-----------------------|--|
| Product Type:         | Primary Antibodies   |
| Clone Name:           | B418M  |
| Applications:         | IF, Neutralize, WB   |
| Recommended Dilution: | <b>ELISA.</b><br><b>Western Blot.</b><br><b>Has Toxin Neutralization activity (<i>in vivo</i> assay).</b>                                      |
| Host:                 | Mouse  |
| Isotype:              | IgG2a  |
| Clonality:            | Monoclonal   |
| Immunogen:            | Supernatant from <i>C. tetani</i>  |
| Specificity:          | This antibody reacts with Tetanus Toxin.   |
| Formulation:          | PBS, pH 7.4<br>State: Purified<br>State: Liquid purified Ig fraction from Ascites (> 95% pure by SDS-PAGE)<br>Preservative: 0.09% Sodium Azide |
| Concentration:        | lot specific   |
| Purification:         | Affinity Chromatography on Protein G   |
| Conjugation:          | Unconjugated   |
| Storage:              | Store the antibody undiluted at 2-8°C.   |
| Stability:            | Shelf life: one year from despatch.  |
| Database Link:        | <a href="#">P04958</a>   |
| Synonyms:             | Clostridium tetani toxin, tetX   |



[View online »](#)