

Product datasheet for **AM31431PU-N**

CTC_RS14060 Mouse Monoclonal Antibody [Clone ID: B416M]

Product data:

Product Type:	Primary Antibodies
Clone Name:	B416M
Applications:	IF, Neutralize, WB
Recommended Dilution:	ELISA. Western Blot.
Reactivity:	Clostridium tetanis
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Supernatant from C. tetani
Specificity:	This antibody reacts with tetanus toxin. Has toxin neutralization activity (in vivo assay).
Formulation:	PBS, pH 7.4 containing 0.09% Sodium azide as preservative State: Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Affinity chromatography on Protein G (>95% pure (SDS-PAGE))
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.
Database Link:	P04958
Synonyms:	Clostridium tetani toxin, tetX



[View online »](#)