

Product datasheet for **AM31434PU-N**

CKMB Mouse Monoclonal Antibody [Clone ID: B927M]

Product data:

Product Type:	Primary Antibodies
Clone Name:	B927M
Applications:	ELISA, LF
Recommended Dilution:	ELISA. Lateral Flow.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Human CK-MB
Specificity:	This antibody recognizes Human Creatine Kinase MB (CK-MB). Other species not tested.
Formulation:	15mM Potassium phosphate, pH 7.3 containing 0.15M Sodium Chloride State: Purified State: Liquid purified Ig fraction Preservative: 0.05% Sodium Azide
Concentration:	lot specific
Purification:	Affinity Chromatography on Protein G
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.
Synonyms:	CK-MB, Creatine Kinase MB



[View online »](#)