

Product datasheet for **AM31437PU-N**

Cyfra-21-1 Mouse Monoclonal Antibody [Clone ID: B974M]

Product data:

Product Type:	Primary Antibodies
Clone Name:	B974M
Applications:	ELISA, WB
Recommended Dilution:	ELISA. Western Blot.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Soluble cytokeratin antigen fraction purified from human cell line
Specificity:	This antibody reacts to Cyfra-21-1, Cytokeratin 19 fragment.
Formulation:	PBS containing 0.05% Sodium azide as preservative State: Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Ion exchange chromatography (>90% pure (SDS-PAGE))
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.
Synonyms:	Cytokeratin 19 fragment



[View online »](#)