

## Product datasheet for **AM31464PU-N**

### FITC Mouse Monoclonal Antibody [Clone ID: 19F/2A3]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	19F/2A3
Applications:	ELISA, IHC
Recommended Dilution:	Detection of FITC-labeled antigens or antibodies. Capture antibody for FITC-labeled biopolymers. Can be enzyme-labeled for use as a secondary antibody for immunohistochemistry with FITC-conjugated mouse monoclonal antibodies.
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	FITC-BSA conjugate
Specificity:	This antibody recognizes Fluorescein and FITC-labeled protein. Reacts with both free and conjugated FITC.
Formulation:	PBS, pH 7.4 containing 0.09% Sodium azide as preservative State: Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Affinity chromatography on Protein A (>90% pure (SDS-PAGE))
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.
Synonyms:	Fluorescein Isothiocyanate



[View online »](#)