

Product datasheet for **AM31465PU-N**

FSH beta (FSHB) Mouse Monoclonal Antibody [Clone ID: 090-17275]

Product data:

Product Type:	Primary Antibodies
Clone Name:	090-17275
Applications:	ELISA
Recommended Dilution:	Suitable for use in ELISA . Recommended antibody pairs for Sandwich Immunoassay: <u>Capture / Detection:</u> AM31465PU-N / BM2500 AM31465PU-N / AM31467PU-N
Reactivity:	Human
Host:	Mouse
Isotype:	IgG
Clonality:	Monoclonal
Specificity:	This antibody recognizes Follicle Stimulating Hormone (FSH).
Formulation:	10mM Phosphate, pH 7.4 containing 150mM Sodium Chloride and 0.09% Sodium Azide as preservative. State: Purified State: Liquid purified Ig fraction (>90% pure by SDS-PAGE). Product is 0.2µm filtered.
Concentration:	lot specific
Purification:	Protein A Chromatography.
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one week or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	follicle stimulating hormone beta subunit
Database Link:	Entrez Gene 2488 Human P01225



[View online »](#)

Background:	<p>FSH is a pituitary hormone involved in the maturation of ovarian follicles and estrogen secretion in females. In the pituitary gland, FSH is produced by gonadotrophs. In males, FSH stimulates the secretion of testosterone.</p> <p>Follicle stimulating hormone enables ovarian folliculogenesis to the antral follicle stage and is essential for Sertoli cell proliferation and maintenance of sperm quality in the testis.</p> <p>Members of the pituitary glycoprotein hormone family, of which FSH is one (see also luteinizing hormone, chorionic gonadotropin, and thyroid stimulating hormone), consist of a shared alpha chain and a beta chain encoded by a separate gene.</p>
Synonyms:	Follitropin beta chain, FSHB, FSH beta
Note:	Source: Ascites
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways:	GnRH signaling pathway, Neuroactive ligand-receptor interaction