

Product datasheet for **AM31470PU-N**

Giardia lamblia Mouse Monoclonal Antibody [Clone ID: BDI361]

Product data:

Product Type:	Primary Antibodies
Clone Name:	BDI361
Applications:	ELISA, IF, WB
Recommended Dilution:	Western Blot. Immunofluorescence. ELISA. Recommended antibody pairs for Sandwich Immunoassay (Capture / Detection): AM31469PU-N / AM31470PU-N AM00845PU-N / AM31470PU-N
Reactivity:	Giardia lamblia
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Intact formalin-fixed cysts
Specificity:	This antibody is specific to the 65kDa <i>Giardia antigen</i> . Detects cysts.
Formulation:	0.01M PBS, pH 7.2 containing 0.09% Sodium azide as preservative This product contains no stabilizing proteins. State: Purified State: Liquid purified Ig fraction (>90% pure (SDS-PAGE))
Concentration:	lot specific
Purification:	Affinity chromatography on Protein A
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.



[View online »](#)