

Product datasheet for **AM31479AF-N**

Hemoglobin Mouse Monoclonal Antibody [Clone ID: 902]

Product data:

Product Type:	Primary Antibodies
Clone Name:	902
Applications:	ELISA
Recommended Dilution:	ELISA. May be used in immunoassay to detect and quantitate human hemoglobin.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Purified Human Hemoglobin
Specificity:	This antibody reacts to Human Hemoglobin. Does not cross-react with hemoglobin of the following species: bovine, chicken, horse, pig and sheep.
Formulation:	PBS, pH 7.4 State: Azide Free State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Affinity chromatography on Protein G
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C to -70°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	Haemoglobin



[View online »](#)