

## Product datasheet for **AM31502PU-N**

### Hepatitis B Core Antigen / HBcAg Mouse Monoclonal Antibody [Clone ID: B945M]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	B945M
Applications:	ELISA
Recommended Dilution:	<b>ELISA.</b> <u>Recommended antibody pair for sandwich immunoassay:</u> Capture Detection AM31501PU-N AM31502PU-N
Reactivity:	Hepatitis B Virus
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Full-length HBcAg recombinant
Specificity:	This antibody reacts to Hepatitis B Virus Core Antigen (HBcAg).
Formulation:	0.01M PBS, pH 7.2 containing 0.09% Sodium azide as preservative This product contains no stabilizing proteins State: Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Affinity chromatography on Protein A (>90% pure (SDS-PAGE))
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	HBV Capsid protein, HBV Core protein, p21.5



[View online »](#)