

Product datasheet for **AM31503PU-N**

Hepatitis B E Antigen / HBeAg Mouse Monoclonal Antibody [Clone ID: BDI831]

Product data:

Product Type:	Primary Antibodies
Clone Name:	BDI831
Applications:	ELISA
Recommended Dilution:	ELISA. <u>Recommended pair for sandwich immunoassay:</u> Capture Detection AM00875PU-N AM31503PU-N
Reactivity:	Hepatitis B Virus
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Specificity:	This antibody reacts to Hepatitis B 'e' Antigen (HBeAg).
Formulation:	0.01M PBS, pH 7.2 containing 0.09% Sodium azide as preservative This product contains no stabilizing proteins. State: Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Affinity chromatography on Protein A (>90% pure (SDS-PAGE))
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)