

## Product datasheet for **AM31509PU-N**

### Hepatitis B Surface Antigen / HBsAg (ad/ay) Mouse Monoclonal Antibody [Clone ID: B521M]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	B521M
Applications:	ELISA, LF
Recommended Dilution:	ELISA. Lateral Flow Assay.
Reactivity:	Hepatitis B Virus
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Purified Hepatitis B Surface Antigen
Specificity:	This antibody reacts to to Hepatitis B Surface Antigen (HBsAg) (ad/ay).
Formulation:	Borate saline buffer containing 0.09% Sodium azide as preservative State: Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Ion exchange chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	HBV surface antigen, Hepatitis B Virus



[View online »](#)