

## Product datasheet for **AM31550PU-N**

### Human IgE Mouse Monoclonal Antibody [Clone ID: 904]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgE Mouse Monoclonal Antibody [Clone ID: 904]
Clone Name:	904
Applications:	ELISA
Recommended Dilution:	ELISA.
Reactivity:	Human
Host:	Mouse
Immunogen:	Human Myeloma IgE
Isotype:	IgG2b
Formulation:	0.015M Potassium phosphate buffer, 0.15M Sodium chloride, pH 7.2 containing 50% glycerol and 0.09% Sodium azide State: Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	DEAE-column Chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)