

Product datasheet for **AM31552PU-N**

Human IgG Mouse Monoclonal Antibody [Clone ID: B573M]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG Mouse Monoclonal Antibody [Clone ID: B573M]
Clone Name:	B573M
Applications:	ELISA, LF
Recommended Dilution:	ELISA. Lateral Flow Assays.
Reactivity:	Human
Host:	Mouse
Immunogen:	Purified Human IgG
Isotype:	IgG1
Formulation:	PBS containing 0.09% Sodium azide as preservative State: Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Ion exchange chromatography (91.02%)
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for 2 weeks or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)