

## Product datasheet for **AM31553PU-N**

### Human IgG (Fc specific) Mouse Monoclonal Antibody [Clone ID: MK1A6]

#### Product data:

|                       |                                                                                                                                             |
|-----------------------|---------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:         | Secondary Antibodies                                                                                                                        |
| Product Name:         | Human IgG (Fc specific) Mouse Monoclonal Antibody [Clone ID: MK1A6]                                                                         |
| Clone Name:           | MK1A6                                                                                                                                       |
| Applications:         | ELISA, IHC                                                                                                                                  |
| Recommended Dilution: | ELISA assay development of total IgG.<br>Passive hemagglutination assays.<br>Histological staining of Acetone Fixed Frozen tissue sections. |
| Reactivity:           | Human                                                                                                                                       |
| Host:                 | Mouse                                                                                                                                       |
| Immunogen:            | Human IgM Rheumatoid Factor                                                                                                                 |
| Isotype:              | IgG1                                                                                                                                        |
| Formulation:          | PBS, pH 7.1<br>State: Purified<br>State: Liquid purified Ig fraction from Ascites<br>Preservative: 0.1mM Sodium Azide                       |
| Concentration:        | lot specific                                                                                                                                |
| Purification:         | Affinity chromatography on Protein A                                                                                                        |
| Conjugation:          | Unconjugated                                                                                                                                |
| Storage:              | Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.<br>Avoid repeated freezing and thawing.                        |



[View online »](#)