

## Product datasheet for **AM31555PU-N**

### Human IgG (Fc specific) Mouse Monoclonal Antibody [Clone ID: 21G/3D3]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG (Fc specific) Mouse Monoclonal Antibody [Clone ID: 21G/3D3]
Clone Name:	21G/3D3
Applications:	ELISA
Recommended Dilution:	<b>ELISA.</b>
Reactivity:	Human
Host:	Mouse
Immunogen:	IgG from human serum
Isotype:	IgG
Formulation:	PBS, pH 7.4 containing 0.09% Sodium azide as preservative State: Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Affinity chromatography on Protein G (>95% pure by SDS-PAGE)
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C.



[View online »](#)