

## Product datasheet for **AM31556AF-N**

### Human IgG (heavy chain) Mouse Monoclonal Antibody [Clone ID: BDI412]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG (heavy chain) Mouse Monoclonal Antibody [Clone ID: BDI412]
Clone Name:	BDI412
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA. Western Blot. Immunohistochemistry on frozen sections.
Reactivity:	Human
Host:	Mouse
Immunogen:	Purified intact human IgG
Isotype:	IgG1
Formulation:	0.15M PBS, pH 7.4 State: Azide Free State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Affinity chromatography on Protein A
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C to -70°C. Avoid repeated freezing and thawing.



[View online »](#)