

Product datasheet for **AM31567PU-N**

Influenza A H5N1 (Haemagglutinin) Mouse Monoclonal Antibody [Clone ID: B253M]

Product data:

Product Type:	Primary Antibodies
Clone Name:	B253M
Applications:	ELISA
Recommended Dilution:	ELISA. <u>Recommended pairs for sandwich immunoassay:</u> Capture Detection AM31567PU-N AM31566PU-N AM31567PU-N AM00942PU-N AM31567PU-N AM00944PU-N
Reactivity:	Influenza A Virus
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Purified avian influenza virus type A (H5N1)
Specificity:	This antibody reacts to Influenza A Hemagglutinin H5 (Avian H5N1).
Formulation:	PBS, pH 7.4 containing 0.09% Sodium azide as preservative State: Aff - Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Affinity chromatography on Protein G (>95% pure (SDS-PAGE))
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.
Synonyms:	Avian Influenza A H5N1 H5 Hemagglutinin



[View online »](#)