

Product datasheet for **AM31579PU-N**

Ketamine Mouse Monoclonal Antibody [Clone ID: B924M]

Product data:

Product Type:	Primary Antibodies
Clone Name:	B924M
Applications:	ELISA, IF
Recommended Dilution:	ELISA. Immunofluorescence.
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Ketamine-BgG
Specificity:	This antibody reacts with Ketamine.
Formulation:	0.05M Phosphate, pH 7.15 containing 0.85% Sodium Chloride State: Purified State: Liquid purified IgG fraction from Ascites Preservative: 0.05% Sodium Azide
Concentration:	lot specific
Purification:	Affinity Chromatography on Protein A
Conjugation:	Unconjugated
Storage:	Store undiluted at 2–8°C for one month or (in aliquots) at –20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.

