

Product datasheet for **AM31589PU-N**

Mycoplasma hominis Surface antigen p120 Mouse Monoclonal Antibody [Clone ID: B565M]

Product data:

Product Type:	Primary Antibodies
Clone Name:	B565M
Applications:	ELISA, IF
Recommended Dilution:	ELISA. Immunofluorescence.
Reactivity:	Mycoplasma hominis
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Recombinant M. hominis 120 kDa Major Surface Protein
Specificity:	This antibody recognizes 120kDa Major Surface Protein of Mycoplasma hominis. Does not cross react with M. pneumoniae and U. urealyticum.
Formulation:	0.01M PBS, pH 7.2 containing 0.09% Sodium azide as preservative This product contains no stabilizing proteins. State: Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Affinity chromatography on Protein G (>95% pure (SDS-PAGE))
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)