

Product datasheet for **AM31594PU-N**

Mycoplasma pneumoniae Mouse Monoclonal Antibody [Clone ID: B749M]

Product data:

Product Type:	Primary Antibodies
Clone Name:	B749M
Applications:	ELISA, IF, WB
Recommended Dilution:	ELISA. Western Blot. Western Blot.
Reactivity:	Mycoplasma pneumoniae
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Partially-purified Mycoplasma pneumoniae
Specificity:	This antibody reacts with multiple bands of Mycoplasma pneumoniae on Western blot. Negative with uninfected cells.
Formulation:	0.01M PBS, pH 7.2 containing 0.09% Sodium azide as preservative This product contains no stabilizing proteins. State: Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Affinity chromatography on Protein A (>90% pure (SDS-PAGE))
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)