

## Product datasheet for **AM31613PU-N**

### Oxycodone Mouse Monoclonal Antibody [Clone ID: B879M]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	B879M
Applications:	ELISA, LF
Recommended Dilution:	Lateral flow. ELISA.
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Oxycodone-BSA Conjugate
Specificity:	This antibody reacts to Oxycodone.
Formulation:	PBS, pH 7.0 containing 0.02% Sodium azide as preservative State: Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	>95% pure (SDS-PAGE)
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)