

Product datasheet for **AM31627PU-N**

Progesterone Mouse Monoclonal Antibody [Clone ID: 254-12113]

Product data:

Product Type:	Primary Antibodies
Clone Name:	254-12113
Applications:	ELISA
Recommended Dilution:	ELISA.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Progesterone-BSA (linked at the 3 position)
Specificity:	This antibody reacts to Progesterone. Cross-reactivity: <i>Test Compound / Progesterone</i> 1µg cortisol =5ng 10ng estradiol <<5ng 500ng estriol <<5ng 20ng testosterone <<5ng
Formulation:	PBS, pH 7.15 containing 0.05% Sodium azide as preservative State: Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Affinity chromatography on Protein A
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one week or (in aliquots) at -20°C for longer. If aliquoted for long term storage, fill volume should be equal to or greater than 50% of the nominal fill volume of the vial used. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	20-dione, pregn-4-ene-3



[View online »](#)