

Product datasheet for **AM31679PU-N**

Tobramycin Mouse Monoclonal Antibody [Clone ID: 128-10052]

Product data:

Product Type:	Primary Antibodies
Clone Name:	128-10052
Applications:	ELISA
Recommended Dilution:	ELISA.
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Tobramycin BSA
Specificity:	This antibody reacts to Tobramycin.
Formulation:	10mM Phosphate, pH 7.4 containing 150mM Sodium chloride and 0.09% Sodium azide State: Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Affinity chromatography on Protein A
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one week or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	Tobrex, Tobracin, Tenebrimycin, Tobramitsetin, Tobramycetin, Gotabiotic, Tobradistin, Tobramaxin, Tobramicin



[View online »](#)