

Product datasheet for **AM31680PU-N**

Toxoplasma gondii p30 Mouse Monoclonal Antibody [Clone ID: B754M]

Product data:

Product Type:	Primary Antibodies
Clone Name:	B754M
Applications:	ELISA, IF, WB
Recommended Dilution:	ELISA. Western Blot. Immunofluorescence.
Reactivity:	Toxoplasma gondii
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Native T. gondii lysate
Specificity:	This antibody reacts to Toxoplasma gondii, SAG1.
Formulation:	0.01M PBS, pH 7.2 containing 0.09% Sodium azide as preservative This product contains no stabilizing proteins. State: Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Affinity chromatography on Protein A (>90% pure (SDS-PAGE))
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Database Link:	Q27298
Synonyms:	SAG1



[View online »](#)