

Product datasheet for **AM31681PU-N**

Tramadol Mouse Monoclonal Antibody [Clone ID: B923M]

Product data:

Product Type:	Primary Antibodies
Clone Name:	B923M
Applications:	ELISA
Recommended Dilution:	ELISA.
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Tramadol-BgG
Specificity:	This antibody reacts to Tramadol. <u>Cross-reactivity:</u> (Based on ELISA(Tramadol at 10ng/ml)) (+/-) Chlorpheniramine 0.1% Diphenhydramine 0.3% Pheniramine 0.04% PCM 0.2%
Formulation:	0.05M Phosphate, pH 7.15 containing 0.85% Sodium chloride and 0.05% Sodium azide as preservative State: Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Affinity chromatography on Protein A
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)