

Product datasheet for **AM31682PU-N**

Trichomonas foetus Mouse Monoclonal Antibody [Clone ID: B1213M]

Product data:

Product Type:	Primary Antibodies
Clone Name:	B1213M
Applications:	ELISA, IF
Recommended Dilution:	Immunofluorescence. ELISA. <u>Recommended antibody pair for sandwich immunoassay:</u> Capture Detection AM31682PU-N AM31683PU-N
Reactivity:	Trichomonas foetus
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Immunogen:	Disrupted Trichomonas foetus
Specificity:	This antibody reacts with <i>Trichomonas foetus</i> . Does not cross react with <i>Trichomonas vaginalis</i> .
Formulation:	0.01M PBS, pH 7.2 containing 0.09% Sodium azide as preservative State: Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Affinity chromatography on Protein A (>90% pure (SDS-PAGE))
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for 2 weeks or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)