

Product datasheet for **AM31688PU-N**

Trichomonas vaginalis Mouse Monoclonal Antibody [Clone ID: B984M]

Product data:

Product Type:	Primary Antibodies
Clone Name:	B984M
Applications:	IF
Recommended Dilution:	Immunofluorescence. Not recommended for use in ELISA.
Reactivity:	Trichomonas vaginalis
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Disrupted Trichomonas vaginalis
Specificity:	This antibody reacts to <i>Trichomonas vaginalis</i> .
Formulation:	0.01M PBS, pH 7.2 containing 0.09% Sodium azide as preservative State: Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Affinity chromatography on Protein A (>90% pure (SDS-PAGE))
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for two weeks or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)