

Product datasheet for **AM31834RP-N**

C3ar1 Mouse Monoclonal Antibody [Clone ID: 74]

Product data:

Product Type:	Primary Antibodies
Clone Name:	74
Applications:	FC, IHC
Recommended Dilution:	Flow Cytometry. Immunohistochemistry on Frozen Sections.
Reactivity:	Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	RBL-2H3 transfectants expressing rat C3aR <u>Donor:</u> BALB/c <u>Fusion Partner:</u> X63-Ag8.653
Specificity:	This Rat C3a receptor antibody detects a recombinant peptide corresponding to amino acids 161-321 of the large extracellular loop structure, which shares 34% homology with Mouse and Human. Deduced amino acid sequences of Human/Rat C3aR is 58.5% and 90.7% with Rat/Mouse.
Formulation:	PBS containing 0.02% Sodium Azide as preservative and EIA grade BSA was added as a stabilizing protein to bring the total protein concentration to 4-5 mg/ml. Label: PE State: Liquid purified Ig fraction Label: R-Phycoerythrin
Concentration:	lot specific
Purification:	Affinity Chromatography on Protein G
Conjugation:	PE
Storage:	Store the antibody undiluted at 2-8°C. DO NOT FREEZE!
Stability:	Shelf life: one year from despatch.
Gene Name:	complement component 3a receptor 1



[View online »](#)

Database Link: [Entrez Gene 84007 Rat O55197](#)

Background: C3aR binds to anaphylotoxin C3a which, among other things, is involved in mast cell degranulation and recruitment of immune cells to the site of inflammation. Activation of the complement cascade produces various fragments, including C3a and C5a. C3aR is characterized by seven transmembrane domains including a large extracellular loop which is coupled to a Gi protein. C3aR expression is detected on glial cells, neurons and cells belonging to the mononuclear phagocyte system. C3aR expression was shown to moderately increase after LPS stimulation.

Synonyms: C3a-R, C3AR, C3AR1, AZ3B, C3R1, HNFAG09

Note:

Protocol: **FLOW CYTOMETRY ANALYSIS:**

Method:

1. Prepare a cell suspension in media A. For cell preparations, deplete the red blood cell population with Lympholyte®-Rat cell separation medium.
2. Wash 2 times.
3. Resuspend the cells to a concentration 2×10^7 cells/ml in media A. Add 50 μ l of this suspension to each tube (each tube will then contain 1×10^6 cells, representing 1 test).
4. To each tube, add 0.5 μ g of this antibody.
5. Vortex the tubes to ensure thorough mixing of antibody and cells.
6. Incubate the tubes for 30 minutes at 4°C.
7. Wash 2 times at 4°C.
8. Add 100 μ l of secondary antibody (FITC Goat anti-mouse IgG) at 1/500 dilution.
9. Incubate the tubes at 4°C for 30-60 minutes.
(It is recommended that the tubes are protected from light since most fluorochromes are light sensitive).
10. Wash 2 times at 4°C in media B.
11. Resuspend the cell pellet in 50 μ l ice cold media B.
12. Transfer to suitable tubes for flow cytometric analysis containing 15 μ l of propidium iodide at 0.5 mg/ml in PBS. This stains dead cells by intercalating in DNA.

Media:

- A. Phosphate buffered saline (pH 7.2) + 5% normal serum of host species + sodium azide (100 μ l of 2M sodium azide in 100 mls).
- B. Phosphate buffered saline (pH 7.2) + 0.5% Bovine serum albumin + sodium azide (100 μ l of 2M sodium azide in 100 mls).

Results:

Tissue Distribution by Flow Cytometry Analysis:

Rat Strain: Wistar

Cell Concentration: 1×10^6 cells per test

Antibody Concentration Used: 2.0 μ g / 10^6 cells

Isotypic Control: FITC Mouse IgG1

Cell Source:

Peritoneal Macrophages 97.2%

Bone Marrow 81.5%