

Product datasheet for **AM31855BT-N**

Human IgG (Fd region) Mouse Monoclonal Antibody [Clone ID: HP6045]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG (Fd region) Mouse Monoclonal Antibody [Clone ID: HP6045]
Clone Name:	HP6045
Applications:	ELISA, R, WB
Recommended Dilution:	ELISA (both capture and detection). Western Blot. Radioimmunoassay
Reactivity:	Human
Host:	Mouse
Immunogen:	Human IgG Fd
Isotype:	IgG1
Formulation:	PBS containing 0.02% Sodium Azide as preservative and EIA grade BSA as a stabilizing protein to bring total protein concentration to 4-5 mg/ml. Label: Biotin State: Liquid purified IgG fraction.
Concentration:	lot specific
Conjugation:	Biotin
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)