

## Product datasheet for **AM31855FC-N**

### Human IgG (Fd region) Mouse Monoclonal Antibody [Clone ID: HP6045]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG (Fd region) Mouse Monoclonal Antibody [Clone ID: HP6045]
Clone Name:	HP6045
Applications:	ELISA, R, WB
Recommended Dilution:	<b>ELISA</b> (both capture and detection). <b>Western Blot.</b> <b>Radioimmunoassay</b>
Reactivity:	Human
Host:	Mouse
Immunogen:	Human IgG Fd
Isotype:	IgG1
Formulation:	PBS containing 0.02% Sodium Azide as preservative and EIA grade BSA as a stabilizing protein to bring total protein concentration to 4-5 mg/ml. Label: FITC State: Liquid purified IgG fraction.
Concentration:	lot specific
Conjugation:	FITC
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)