

## Product datasheet for **AM31855HR-N**

### Human IgG (Fd region) Mouse Monoclonal Antibody [Clone ID: HP6045]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG (Fd region) Mouse Monoclonal Antibody [Clone ID: HP6045]
Clone Name:	HP6045
Applications:	ELISA, R, WB
Recommended Dilution:	<b>ELISA</b> (both capture and detection). <b>Western Blot.</b> <b>Radioimmunoassay.</b>
Reactivity:	Human
Host:	Mouse
Immunogen:	Human IgG Fd
Isotype:	IgG1
Formulation:	HRP conjugated IgG buffered in PBS with 40% Glycerol and EIA grade BSA as a stabilizing protein to bring the total protein concentration to 4-5 mg/ml. No Sodium Azide with other preservatives. Label: HRP State: Liquid purified IgG fraction Preservative: None
Concentration:	lot specific
Conjugation:	HRP
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)