

## Product datasheet for AM31856BT-N

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **BCL3 Hamster Monoclonal Antibody [Clone ID: Ham150-3.5]**

**Product data:** 

**Product Type:** Primary Antibodies

Clone Name: Ham150-3.5

**Applications:** ELISA

**Recommended Dilution:** This Clone has been described to work in**ELISA**, **Western Blot** and **Flow Cytometry**.

Reactivity: Human, Mouse

Host: Hamster Isotype: IgG1

Clonality: Monoclonal

**Immunogen:** Recombinant Bcl-3-6-histidine tag from Armenian Hamster.

Fusion Partner: SP2/0 myeloma.

**Specificity:** Recognizes B-cell Leukemia/lymphoma 3

Formulation: PBS containing 0.02% Sodium Azide as preservative and EIA grade BSA as a stabilizing protein

to bring total protein concentration to 4-5 mg/ml.

Label: Biotin

State: Liquid purified IgG fraction.

**Concentration:** lot specific

**Purification:** Protein G Affinity Chromatography.

**Conjugation:** Biotin

Storage: Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.

Avoid repeated freezing and thawing.

**Stability:** Shelf life: one year from despatch.

Gene Name: B-cell CLL/lymphoma 3

Database Link: Entrez Gene 602 Human

P20749





## BCL3 Hamster Monoclonal Antibody [Clone ID: Ham150-3.5] - AM31856BT-N

**Background:** B cell Lymphoma 3 is a member of the IkB subfamily of inhibitors. This subfamily is part of the

NF-kB transcription factor protein family and it is suggested that a balance in the

concentration of various NF-kB family members regulates apoptosis and survival of activated

T cells.

Bcl-3 has been shown to have versatile functions such as cytoplasmic activation of p50 homodimers, translocation to the nucleus, and modulation of the transcriptional machinery in the nucleus. Bcl-3 activity often relies on several nuclear interacting proteins such as Tip60, Jab1, Bard1, and Pirin. The Mouse Bcl-3 coding region exhibits 80% homology with Human Bcl-3 which is associated with human B-cell chronic lymphocytic leukemias (CLLs).

Bcl-3 is detected in various tissues such as spleen and other lymphoid organs.

Synonyms: BCL3, BCL4, D19S37

**Protein Families:** Druggable Genome, Transcription Factors