

## Product datasheet for **AM31863RP-N**

### TCR DO-11.10 Mouse Monoclonal Antibody [Clone ID: KJ1-26]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	KJ1-26
Applications:	FC
Recommended Dilution:	<b>Flow Cytometry..</b> This clone has also been reported to work in <b>Immunoprecipitation, ELISA</b> and <b>Inhibition Assay.</b>
Reactivity:	Mouse
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	T cell hybridoma, D011.10. Donor: Mouse ascites fluid.
Specificity:	This Monoclonal Antibody binds the T-cell receptor (TCR) expressed on DO11.10 and DO11.10.24 murine T cell hybridomas, but not on any other T cell hybridomas tested. (1) The TCR expressed on DO11.10 recognizes an ovalbumin peptide in the context of I-Ad major histocompatibility molecules. (1,2) The use of DO11.10 and Clone KJ1-26 have proven to be instrumental in describing many of the intricate details of T cell immunity. (3-6) The DO11.10 T cell hybridoma is available from the laboratories of Drs. Kappler and Marrack. (2,7)
Formulation:	PBS containing 0.02% Sodium Azide as preservative and EIA grade BSA as a stabilizing protein to bring total protein concentration to 4-5 mg/ml. Label: PE State: Liquid purified IgG fraction.
Concentration:	lot specific
Conjugation:	PE
Storage:	Store the antibody undiluted at 2-8°C. <b>DO NOT FREEZE!</b> Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)

**Synonyms:** DO-11.10 T cell Receptor