

Product datasheet for **AM31884BT-L**

CD11a / ITGAL Rat Monoclonal Antibody [Clone ID: IBL-6/2]

Product data:

Product Type:	Primary Antibodies
Clone Name:	IBL-6/2
Applications:	FC
Recommended Dilution:	Flow Cytometry.
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	EL4 (mouse thymoma cell line) <u>Donor:</u> Wistar spleen <u>Fusion Partner:</u> Sp-2/0Ag14
Specificity:	This antibody reacts to CD11a / ITGAL.
Formulation:	Label: Biotin State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Affinity chromatography on Protein G
Conjugation:	Biotin
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Database Link:	<u>P24063</u>
Background:	The CD11a is expressed on leukocytes and it mediates cell-cell and cell-adhesion by binding to CD54 and ICAM-2).
Synonyms:	Integrin alpha-L, LFA1, LFA-1



[View online »](#)